

# Schedule C Planned Maintenance

Site: <u>Green Light Imaging</u>	Month: <u>May 2021</u>
System ID#: <u>Mobil 12x</u>	Dispatch/Job #: _____
Service Engineer: <u>Michael Stephens</u>	

PM Check	Test Purpose	Pass	Fail	Comments
<b>Image Quality</b>				
<u>Alignment Light Check</u>	IQ	✓		
<u>LVShim</u>	IQ	✓		Adjust
<u>Grafidy3</u>	IQ			N/A
<u>Coherent Noise</u>	IQ	✓		
<u>Signal to Noise Check</u>	IQ	✓		
<u>Spike Noise Check</u>	IQ	✓		
<u>PM Check (Service Contract Customers Only)</u>	IQ	✓		

PM Check	Test Purpose	Pass	Fail	Comments
<b>Magnet Room</b>				
<u>Inspect Front cover Cable Take-up (Saltbox style magnets only)</u>	check			N/A
<u>PAC Leakage Current Test</u>	Safety	✓		
<u>Cardiac Gating Cable</u>	Safety	✓		
<u>Patient Blower &amp; Filter</u>	Check	✓		cleaned
<b>RF/System Cabinet</b>				
<u>Check/clean filters in RFS cabinet</u>	Clean	✓		
<u>Fans &amp; Filter for SRFD</u>	Clean	✓		
<u>Check RF Output Power</u>	Check	✓		
<b>Patient Handling</b>				
<u>Patient Table Checks</u>	Safety	✓		Adj
<b>Gradients</b>				
<u>Fans &amp; Filter</u>	Clean	✓		
<u>Clean Reservoir, condenser fins, &amp; pump strainer (Neslab only)</u>	Clean			N/A
<u>Check Fluid levels and Valve of Heat Exchanger</u>	Check	✓		
<b>Twin Accessory Cabinet</b>				
<u>Filter Replacement</u>	Replace			N/A
<u>Clean Inlet Screen</u>	Clean			
<u>Solenoid Valve Replacement</u>	Replace			
<b>Computer</b>				
<u>Storelog</u>	Check	✓		cleaned
<b>Brainwave Cabinet</b>				
<u>Clean filters on Brainwave Lite PC</u>	Clean			N/A
<b>Magnet</b>				
<u>Test GE Magnet Rundown Unit (MRU)</u>	Safety	✓		pr

Schedule C Signa Planned Maintenance.

Date Completed: 5/28/2021